

TSC22D3-TSC22 domain family, member 3 Gene

Reference: PTXBC018148

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DIP; DSIPI; TSC-22R; TSC22 domain family protein 3; DSIP-immunoreactive leucine zipper protein; DSIP-immunoreactive peptide; TSC-22 related protein; TSC-22-like protein; delta sleep-inducing peptide immunoreactor; glucocorticoid-induced leucine zipper pro		

For research use only.

GENE DESCRIPTION	TSC22 domain family, member 3
GENE ID	1831
ORIGIN SPECIES	Human
ACCESSIONS	BC018148
VECTOR	pDONR223
NCBI SYMBOL	TSC22D3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatctggtagaatacatctgatgtatgtgagagagaggtggagatcctgaaggagcagatccgagagctggaggagaagaactccagctagagcgtgagaacacccgtgtgaagaccctggcaagcccagagcagctggagaagtccagtcctgtgtgagccctgaagagccagctcccgaatcccacaagtgcccaggccctgggtggttctgcggtgttaa