

NDUFB3-NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 3, 12kDa Gene

Reference: PTXBC018183

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B12; CI-B12; NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 3; NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 3, 12kDa; NADH-ubiquinone oxidoreductase B12 subunit; complex I-B12; NADH:ubiquinone oxidoreductase subunit B3		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 3, 12kDa
GENE ID	4709
ORIGIN SPECIES	Human
ACCESSIONS	BC018183
VECTOR	pDONR223
NCBI SYMBOL	NDUFB3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccatgaacatggacatgacatgacataaaatggaactccagattatagacaatggaatagaaggacaccattagaaactccagaagaagctggctgcaaaaggcctaaggatccatggggccgcaatgaagctggagatacatgggtggcttgcgaagagt
gtttccctttctgatgtattcttaaaaggattcaaatggggatttgcctgatttggtagctgtaggagctgaatattacctggagtccttgaataaagataaagaagcatcactga