

MT1X-metallothionein 1X Gene

Reference: PTXBC018190

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MT-11; MT1; metallothionein-1X; MT-1X; MT-IX; metallothionein-IX; testicular tissue protein Li 124; metallothionein 1X		

For research use only.

GENE DESCRIPTION	metallothionein 1X
GENE ID	4501
ORIGIN SPECIES	Human
ACCESSIONS	BC018190
VECTOR	pDONR223
NCBI SYMBOL	MT1X
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacccaactgctcctgctgctgctgtggctcctgtgctgtgccggctcctgcaaatgcaacagtgcaaatgcacctcctgcaagaagagctgtgctcctgctgcccgtgtgggtgtgcccaagtgtgccagggtgcatctgcaaggacgtcagacaagtcagctgctgtg
cctga