

TBCA-tubulin folding cofactor A Gene

Reference: PTXBC018210

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	tubulin-specific chaperone A; CFA; TCP1-chaperonin cofactor A; chaperonin cofactor A; tubulin folding cofactor A		

For research use only.

GENE DESCRIPTION	tubulin folding cofactor A
GENE ID	6902
ORIGIN SPECIES	Human
ACCESSIONS	BC018210
VECTOR	pDONR223
NCBI SYMBOL	TBCA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgatcctcgcgtgagacagatcaagaccggcgtggtgaagcgttggtcaagaaaaagtgtatgataaaaggcgaacaagaagaaagattgaaaaatgagagctgaagacggtgaaaattatgacattaaaaagcaggcagagatcctacaagaatccaggatgatgatccagattgccagcgcagggtggaagcgcataattggtatctcaacggatactagaaaaatgaaaaagacttgaagaagctgaggaatataaagaagcacgtttagtactggattcagtgaaagttagaagcctga