



Reference: PTXBC018269

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KIAA1967; DBC-1; DBC1; NET35; p30 DBC; p30DBC; cell cycle and apoptosis regulator protein 2; cell division cycle and apoptosis regulator protein 2; p30 DBC protein; cell cycle and apoptosis regulator 2		

For research use only.

GENE DESCRIPTION	KIAA1967
GENE ID	57805
ORIGIN SPECIES	Human
ACCESSIONS	BC018269
VECTOR	pDONR223
NCBI SYMBOL	KIAA1967
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]