

ZFP36L1-zinc finger protein 36, C3H type-like 1 Gene

Reference: PTXBC018340

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	mRNA decay activator protein ZFP36L1; BRF1; Berg36; ERF-1; ERF1; RNF162B; TIS11B; cMG1; EGF-response factor 1; TPA-induced sequence 11b; ZFP36-like 1; butyrate response factor 1; early response factor Berg36; zinc finger protein 36, C3H type-like 1; zinc		

For research use only.

GENE DESCRIPTION	zinc finger protein 36, C3H type-like 1
GENE ID	677
ORIGIN SPECIES	Human
ACCESSIONS	BC018340
VECTOR	pDONR223
NCBI SYMBOL	ZFP36L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccaccacctgtgtgtccaccatcttgacttgagcgaagtgttatgaaggtaacaagatgtcaactatagtgctccagtgacgggggttgctgtgtgacagaaaggcagtgaggacccctgtgtgtgggggttccctcggaggcactcagtcaccctgccagctcca
agttccaccagaaccagctcctcagcagcctcaagggtgagccagccccgctctgagctcgcgagacagccgctccgagaccgctccttctcggaagggggcgagcggctgtgtgtccaccagaagcagccggggcgccaggtcaactccagcgtacaagacgggagc
tgtgccgcccccttgaggaaacggtgctgtgtaagtaggggacaagtgccagttcgacacggcatccagagctccgagcctgacccgccaccaccaagtaagacggagctgtgtccgacattccacaccatcggtttgtccctacgggccccgctgccactcatccacaac
gctgaagagcgcctgtgcccgtggccggggccgggacgtctccgctgaccgtccccgctccagcatagcttttagcttgggttccagtgccgctgcacccgctgtgccaccgggctgtgtgacagcccccagtcacacccccccctattctgagcgccgatgacctctg
ggctcactacacctgcccgatggcaccataaccctttgctctccagccaggagctgcaagcctctttgccccatgagatggggtgctcccgaccacctctcttccggccatgtccgagtcacctcacatgtttgactctccccagccctcaggattctctct
cggaccaggaggctacctgagcagctccagcagccacagtggtcagactccccgacctggacaactcaagacgcctgccatctcagcagactttccatctcagatgactaa