



Reference: PTXBC018359

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MAT; MATA1; SAMS; SAMS1; S-adenosylmethionine synthase isoform type-1; MAT 1; MAT-I/III; S-adenosylmethionine synthetase isoform type-1; adoMet synthase 1; adoMet synthetase 1; methionine adenosyltransferase 1; methionine adenosyltransferase I, alpha; met		

For research use only.

GENE DESCRIPTION	methionine adenosyltransferase I, alpha
GENE ID	4143
ORIGIN SPECIES	Human
ACCESSIONS	BC018359
VECTOR	pDONR223
NCBI SYMBOL	MAT1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]