

FGF21-fibroblast growth factor 21 Gene

Reference: PTXBC018404

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	fibroblast growth factor 21		
For research use only.			

GENE DESCRIPTION	fibroblast growth factor 21
GENE ID	26291
ORIGIN SPECIES	Human
ACCESSIONS	BC018404
VECTOR	pDONR223
NCBI SYMBOL	FGF21
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggactcggagacacgggttcgagcactcagggtgtgggtttctgtgctggtcttctgtgggagcctgccaggcacacccatccctgaactccagtcctctcgaattcggggccaagtccggcagcggtacacctacacagatgatgccagcagacagaagccac
ctggagatcaggagagatgggacggtggggggcgctgctgaccagagccccgaaagtctctgcagctgaagccctgaagccgggagttatcaaatcttgggagtcagacatccagggttcctgtgccagcggccagatggggccctgtatgatcgctccacttgacctgaggc
ctgcagcttcgggagctgcttcttgaggacggataaatgtttaccagtcgaaagccacggcctccgctgcacctgccagggaacaagtccccacaccgggacacctgcacccgaggaccagctcgtctcctgccactaccaggcctgccccccgacccccggagccacccgg
aatcctggccccccagccccgatgtgggtctctcgaccctctgagcatgtgggaccttccaggccgaagccccagctacgtctctga