

CRB3-crumbs homolog 3 (Drosophila) Gene

Reference: PTXBC018409

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein crumbs homolog 3; crumbs family member 3; crumbs 3, cell polarity complex component		

For research use only.

GENE DESCRIPTION	crumbs homolog 3 (Drosophila)
GENE ID	92359
ORIGIN SPECIES	Human
ACCESSIONS	BC018409
VECTOR	pDONR223
NCBI SYMBOL	CRB3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgaacccgggctgggctgttctggcgtggcctgccgttctgctggccgctggggccgagcctggggcgaatacagaccacttctgcaaatgagaatagcactgtttgcttccaccagctccagctccgatggcaacctgcgtccggaagccatcactgctatca
tcgtgtgtcttccctcttggctgcttgccttggcgtgtgggctggcactgttgggtgcggaagcttgggagaagcggcagacggagggcacctaccggccagtagcgaggagcagttctccatgcagccgagggccgggcccctcaggactccaaggagacggtgcagggct
gcctgcccctatag