

ZFAND2B-zinc finger, AN1-type domain 2B Gene

Reference: PTXBC018415

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AIRAPL; AN1-type zinc finger protein 2B; AIRAP-like protein; arsenite-inducible RNA-associated protein-like protein; zinc finger, AN1-type domain 2B; zinc finger AN1-type containing 2B		

For research use only.

GENE DESCRIPTION	zinc finger, AN1-type domain 2B
GENE ID	130617
ORIGIN SPECIES	Human
ACCESSIONS	BC018415
VECTOR	pDONR223
NCBI SYMBOL	ZFAND2B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagtttcggacctcggcgtcactgttcggagccgagctgtcagcgttggtttctgccgttaagtgtgatgcctgtcaggcatcttctgcgcagaccatgtggcctacgccagcatcactgtggatctgcttaccaaaaggatatccaggtacctgtgtgccctctctgtaattgtg
cctgtgctgtggccagaggggagccccctgacctgtctgtgggagagcacattgacagagactgtcgtctgtatccagcacagcaaaacgttaagatcttaccataaagtgtgaacgcgtggctgcccggcagcgagaaatgatgaaactgacctgtgaacgtgtagccgaaactt
ctgcatcaagcaccggcatccactggacatgattgtctgtgggagggggcaccacaccgcccggcaggactgtctgcatctccagagcacaagctgtgtgcttctacaagcactgtccccagcccaagtcacacatgcttctctgacctctccagcagagccacaacccgatctcc
gtcctggacagccccctcagtgattgcttgcagaatggcctgagtgaggatgaagctctgcagcgggcccctggaaatgtccctggcagaaacaaacccaggttccaagttgtcaggaggaagaagacctagctttgacaaagcactgtcagccagttaggcagaataaccagcggc
agcagggccagagccgagctcgaagccgtccaactgcagcctgtgctag