



Reference: PTXBC018537

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	solute carrier family 35 member F5; HCV NS5A-transactivated protein 3; hepatitis C virus NS5A-transactivated protein 3		

For research use only.

GENE DESCRIPTION	solute carrier family 35, member F5
GENE ID	80255
ORIGIN SPECIES	Human
ACCESSIONS	BC018537
VECTOR	pDONR223
NCBI SYMBOL	SLC35F5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]