

MAP1LC3B-microtubule-associated protein 1 light chain 3 beta Gene

Reference: PTXBC018634

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MAP1LC3B-a; ATG8F; LC3B; MAP1A/1BLC3; microtubule-associated proteins 1A/1B light chain 3B; MAP1 light chain 3-like protein 2; MAP1A/MAP1B LC3 B; MAP1A/MAP1B light chain 3 B; autophagy-related ubiquitin-like modifier LC3 B; microtubule associated protein		

For research use only.

GENE DESCRIPTION	microtubule-associated protein 1 light chain 3 beta
GENE ID	81631
ORIGIN SPECIES	Human
ACCESSIONS	BC018634
VECTOR	pDONR223
NCBI SYMBOL	MAP1LC3B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgtcgagagaagacctcaagcagcgcgcacctcgaacaagagtagaagatgtccgactattcgagagcagcatccaacaaaatcccggtgataatagaacgatacaagggtgagaagcagcttcctgtctggataaaacaaagttcctgtacctgaccatgtcaacatga
gtgagctcatcaagataaattagaagcgcttacagctcaatgctaatacagcccttcctctgtgtggaacggacacagcatggtcagcgtctccaccaaatactcagaggtgtatgagagtgagaaagatgaagatggattcctgtacatggtctatgcctccaggagacgttcggatgaa
attgtcagtgtaa