



Reference: PTXBC018650

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD363; CHEDG1; D1S3362; ECGF1; EDG-1; EDG1; S1P1; sphingosine 1-phosphate receptor 1; S1P receptor 1; S1P receptor Edg-1; endothelial differentiation G-protein coupled receptor 1; endothelial differentiation, sphingolipid G-protein-coupled receptor, 1; sp		

For research use only.

GENE DESCRIPTION	sphingosine-1-phosphate receptor 1
GENE ID	1901
ORIGIN SPECIES	Human
ACCESSIONS	BC018650
VECTOR	pDONR223
NCBI SYMBOL	S1PR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]