

C16orf59-chromosome 16 open reading frame 59 Gene

Reference: PTXBC018719

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C16orf59; chromosome 16 open reading frame 59		
For research use only.			

GENE DESCRIPTION	chromosome 16 open reading frame 59
GENE ID	80178
ORIGIN SPECIES	Human
ACCESSIONS	BC018719
VECTOR	pDONR223
NCBI SYMBOL	C16orf59
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagccgaacccccaggcctggggcgggctcaggaccagcaaatggccccatccgtctctcaggccccagaagccttcacactcaaggagaaggggacactgtcgggtgctcggcattcaggaaagcagctccagaactcgagcctgtggccagctca
gttccacagaccagtgattccacggatgccgcgtgccaaaaccagttctccagaacatgcagacagcttcaggcggggccccagcccaggtcagtgctgtggagggtggagcggagggggcgcctgcggaaggcctgctcgtgtgagactgcgcatgaggagga
gctctcggcagccccatggactggatgcaggagtagcctgctcctcagctggaggcctgcaggccatgtggggcagtgctgcacagctgcaggagctgcgtgcagcgggtggcggaacagccaccaagaccatgtctgtgggagggccccccggagcctcgcctcc
tgtgggggtagagcggagcctgcatggagccccagctgctgtctactcagcaccagggagctgcagaccctggcggccctcaagctgcgagtggtgtgtgtgacagcagatccacttgaaaaggctcgtgatggtgaactcctccccctggtgaagcgtgcacagccgagg
ggccgcccctggtggccctgtgcccgtgtgcacagcctgctcgcagggaggagcagctgtcttaccatcctgcgggatgaacctgcagctga