

SURF4-surfeit 4 Gene

Reference: PTXBC018741

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ERV29; surfeit locus protein 4; surface 4 integral membrane protein; surfeit 4		

For research use only.

GENE DESCRIPTION	surfeit 4
GENE ID	6836
ORIGIN SPECIES	Human
ACCESSIONS	BC018741
VECTOR	pDONR223
NCBI SYMBOL	SURF4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccagaaacgacctgatggcagggcggaggaactgcgcacagttcctcgtgtcacaaagcagtaacctgccccacgtggcgccctctgtctgacagcaccttctggaggacggcatccgtatgtgttcagtgagcgagcagcgcgactacatcgacaccacctgga
actgcggctacctgctgctcctcctcctcctcaacttgcgggacagctgactggctgcctcctgtgttgagcaggaaacttcgtgagtagcctgctcggcctcttggaaatcatagctctgcagacgattgcctacagcattttagggacttgaagtgttgatgaggaaacctggc
cctgggaggaggcctgttctgctcctagcagaatccgttctgaagggaagagcatgttgcggcgctccccacatgcgtgagagctccccaaacagtagatgcagctcggaggcagggtcttgcgtgtctgatgttcatgacctcttcactttgacgccagctctttctattgtcc
agaacatcgtgggcacagctctgatgatttttagtggccaltggttttaaaaccaagctggcgtttagactctgttgtgtggctcttggcacaacgtatatattcaacgccttctggaccattccagctctacaagccatgcagacttctgaaatacgacttctccagacatgtcgtgattgg
ggccttgcctcctggtggcctgggcccctgggggtgtctccatggatgagaagaagaaggagtggttaa