



Reference: PTXBC018926

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	VRL; VRL-1; VRL1; transient receptor potential cation channel subfamily V member 2; OTRPC2; osm-9-like TRP channel 2; vanilloid receptor-like protein 1		

For research use only.

GENE DESCRIPTION	transient receptor potential cation channel, subfamily V, member 2
GENE ID	51393
ORIGIN SPECIES	Human
ACCESSIONS	BC018926
VECTOR	pDONR223
NCBI SYMBOL	TRPV2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]