



Reference: PTXBC018927

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TSD; beta-hexosaminidase subunit alpha; N-acetyl-beta-glucosaminidase subunit alpha; beta-N-acetylhexosaminidase subunit alpha; hexosaminidase A (alpha polypeptide); hexosaminidase subunit A; hexosaminidase subunit alpha		

For research use only.

GENE DESCRIPTION	hexosaminidase A (alpha polypeptide)
GENE ID	3073
ORIGIN SPECIES	Human
ACCESSIONS	BC018927
VECTOR	pDONR223
NCBI SYMBOL	HEXA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]