



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

PHLDA1-pleckstrin homology-like domain, family A, member 1 Gene

Reference: PTXBC018929

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DT1P1B11; PHRIP; TDAG51; pleckstrin homology-like domain family A member 1; PQ-rich protein; PQR protein; T-cell death-associated gene 51 protein; apoptosis-associated nuclear protein; proline- and glutamine-rich protein; proline- and histidine-rich prote		

For research use only.

GENE DESCRIPTION	pleckstrin homology-like domain, family A, member 1
GENE ID	22822
ORIGIN SPECIES	Human
ACCESSIONS	BC018929
VECTOR	pDONR223
NCBI SYMBOL	PHLDA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr