



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ZNF688-zinc finger protein 688 Gene

Reference: PTXBC018997

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

zinc finger protein 688

For research use only.

GENE DESCRIPTION	zinc finger protein 688
GENE ID	146542
ORIGIN SPECIES	Human
ACCESSIONS	BC018997
VECTOR	pDONR223
NCBI SYMBOL	ZNF688
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagccgcccgcacccgctctggccgagagccggggagaccggcctggfctgaggaagcccgggactfagactctcgagagctggccgctfctactctcccggaggagtggggctgctcgcccccgcgagaggctctgtaccgggcagtgatgcaggagacc
tacggccacctggggcgctcggatccagggccccaacacgccctcatcttggatggaacaggagagtgaggctggagccccgcgccaggaatctgagaaaggggaaagactgggagagagctcgaggagagatgtcccaaacaggaaagagggaaccggagga
agtcccaagagcgaaggcgctagaaagctctctgtaggagagctctgaagctgctggagacgaacctgaccagcatattagctgctggcccccgcacggcagacaccccaaaatgctgctgagaccgaccacagagacacagccccggcacccataaccagatg
gagctcagggccggcagggcgccagctgtgacagctggcgccgcgctctaccctctactctgtgctgacacacagggcagctgacatcggggagcgccctctccccgagctggcagctgcgtctcaagaggaattgcgagtggaagcgccacatggtgatcc
ccctcctctctcgggggggccggagggcgagagcctggatccgctgtgctctggccccctgcagctgacggcgaccgctgctgctctccgacacccagacacacacatctcgagagtgctggcctgga