



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

PRUNE2-prune homolog 2 (Drosophila) Gene

Reference: PTXBC019095

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	truncated PRUNE2; BMCC1; BNIPXL; C9orf65; KIAA0367; protein prune homolog 2; BCH motif-containing molecule at the carboxyl terminal region 1; BNIP2 motif containing molecule at the carboxyl terminal region 1; BNIP2 motif-containing molecule at the C-termi		

For research use only.

GENE DESCRIPTION	prune homolog 2 (Drosophila)
GENE ID	158471
ORIGIN SPECIES	Human
ACCESSIONS	BC019095
VECTOR	pDONR223
NCBI SYMBOL	PRUNE2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagaatttttgcacgcgcgaatactaaactgaatcgaagcaaacgcttgagaaggctcatgtggtttggcgcttaactctgtgacttgattctctacattctacattcacaatagcttactttctagacaaggctcagtcaccagggggttctgtgttaccagtctgcaacataccaaga
actgaattcaactacttcaccgagacgaggtttattttgaagacgttaaattttccgaatcatcacaatatccgggatgaattaaacctgcatacgtataatgatgaagggaagtattcgataaaccttgttgccagcagctgtctgccggagtgaagacaaaactttagaatcagcagttgtca
aagtcattaaaccggttgagacagacgatlcccaactgttgattccgagagcttctctctctctgtctgaagagatcttccaagaggctcctgagctcatcaccgagacaacttgctcatcgctcagaggtagcattctttcaagtggatgaccatggaatcagagaagaatctcagagaa
gcaggaggaaaattcttctatcttggaagaaaaattcttaacttgctccaagagaggacatcatcaactgctcatcaggagaccgattgactgctcaggttgaatgatgaacagacaatgtgaaagatctaaaggagctgcatagggagaataaaagtgccattagtactgtgagc
atgaaccttcaggctaaaggctgacgacgctgctgttttag

More info: Please contact us for any question: contact@proteogenix.fr