

SLC23A1-solute carrier family 23 (nucleobase transporters), member 1 Gene

Reference: PTXBC019225

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SLC23A2; SVCT1; YSPL3; solute carrier family 23 member 1; Na(+)/L-ascorbic acid transporter 1; hSVCT1; sodium-dependent vitamin C transporter-1; solute carrier family 23 (ascorbic acid transporter), member 1; solute carrier family 23 (nucleobase transport		

For research use only.

GENE DESCRIPTION	solute carrier family 23 (nucleobase transporters), member 1
GENE ID	9963
ORIGIN SPECIES	Human
ACCESSIONS	BC019225
VECTOR	pDONR223
NCBI SYMBOL	SLC23A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggggccaggagacgtcgaggggcgacacagcatgaaccaccaggggacccctcgaccccgctaccacagagcctaagttgacatgtgtacaagatcaggagcgtgccaccttggtacctgtgcatcctgtgggcttcagcactacgtgacatgcttcagtggtaccatc
gccgtgcccttctctgctggcgtgagcgctgtgtgtggggccacgaccagcatggttagtcagctcctcggcaccatcttcacgtgcgtgggcatgttctgcactctcttggcatgattacagctgtggggctgtccaacctgcaatttgggacatgaactcctctcgcaacctctctgctg
ggatttccatgttcttcgggctcacgctgcccattacgtgagtcaccaacctggcgccatcaatacaggcattctgaagtggatcagattctgattgtgctgctgaccacggagatgttggggcgggtgcttgccttcatacttgacaacacagtgccaggagagccagagagcgtgg
tctgatacagtggaagctggggctcatgccaacagtgacatgttccagcctcaagagctacgatttcccattgggatgggcatagtaaaagaattaccttctgaaatacttctatctgccagcttcaaaggattttcaagttcaaaagatcagattgcaattccagaagacac
cagaaaatacagaactgcattctgtgtgaccaaggtctga