

PIH1D2-PIH1 domain containing 2 Gene

Reference: PTXBC019238

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PIH1 domain-containing protein 2; PIH1 domain containing 2		
For research use only.			

GENE DESCRIPTION	PIH1 domain containing 2
GENE ID	120379
ORIGIN SPECIES	Human
ACCESSIONS	BC019238
VECTOR	pDONR223
NCBI SYMBOL	PIH1D2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagacatctcaaaaggctgtcttaccgaagtactcagtttggaaacctctgatgatctagctcagagtgacccctgagggtatgagaagttattcagcagcagctgaaagaagggaacagctctgtgctgccccagaaccacagctttgtctacagaccaggattctgaaacaa
aagaaaaatacttttatcaacctgtgtcagtggaagaaggatccagctcccaatcaaccactccagctaccttaactgttggcaaccagaagatacaactgagatatcagatgcttacacagtcattgatgttgcctacaatcctgatgttctcatgcagcagaaaaaggaccaagtga
aaaaaatcagtttaattcagatggccatgaaatgcattgaggagaattccagttcacctctcacactcttaccatattaccaaattagaataaaagggaagcattcaagaatgaaacaaatctgatgggaatccaaactgattccatagatttaagagaaaaaatgagaagggaactaactc
ttggacagatacgaagcagtgactatgagcaatccagatcattctcaactgttactgccaagaccaagtttcaggcgaagcagtggtgtctgatagaagagatttccagtgactgaaatccaggtggagatgaagatccagcctatgaactaaaaattgtgcatgacacagtgagaacc
tctgaaaattgagttgaaagtgaattacctggtatttaattctgtctctctgtgaccttagtgttctgaggatgatttattgattgaagtatctgagaagtacagattacatctgaattctccaaactattgatactgaaatgaccacagcaaaattatcaagaaaaatccacgctaatacaca
atgcctttgtgtga