



Reference: PTXBC019246

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	JM24; SNAT5; pp7194; sodium-coupled neutral amino acid transporter 5; solute carrier family 38 (amino acid transporter), member 5; system N transporter 2; transport system N, protein 2; solute carrier family 38 member 5		

For research use only.

GENE DESCRIPTION	solute carrier family 38, member 5
GENE ID	92745
ORIGIN SPECIES	Human
ACCESSIONS	BC019246
VECTOR	pDONR223
NCBI SYMBOL	SLC38A5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]