

HYPK-Huntingtin interacting protein K Gene

Reference: PTXBC019262

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C15orf63; HSPC136; huntingtin-interacting protein K; huntingtin yeast partner K; huntingtin interacting protein K		

For research use only.

GENE DESCRIPTION	Huntingtin interacting protein K
GENE ID	25764
ORIGIN SPECIES	Human
ACCESSIONS	BC019262
VECTOR	pDONR223
NCBI SYMBOL	HYPK
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccgctgtgaaatagatatggcgaccgaggggatgtggagctggagttggagactgagaccagtggaccagagcggcctccggagaagccacggaaacatgacagcgggtgcggcgacttgagcgggtcaccgactatcagaggagaaggagatccagagttc
caatctggagacggccatgtctgtattggagacagaaggtccgggagcagaaagccaaacaggagcgggagaaagaactggcaaaagtctactatcaagaaggagatctggagctaataatgactgagatggagatatctcgagcagcagcaaacgagtttgcgggaacaca
tgggcaacgtggtagaggcgttattgccctaaccaactga