



Reference: PTXBC019266

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	D19S593E; DMR-N9; DMRN9; gene59; dystrophia myotonica WD repeat-containing protein; dystrophia myotonica-containing WD repeat motif protein; protein 59; dystrophia myotonica, WD repeat containing		

For research use only.

GENE DESCRIPTION	dystrophia myotonica, WD repeat containing
GENE ID	1762
ORIGIN SPECIES	Human
ACCESSIONS	BC019266
VECTOR	pDONR223
NCBI SYMBOL	DMWD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]