

BCL2L1-BCL2-like 1 Gene

Reference: PTXBC019307

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BCL-XL/S; BCL2L; BCLX; Bcl-X; PPP1R52; bcl-2-like protein 1; apoptosis regulator Bcl-X; protein phosphatase 1, regulatory subunit 52; BCL2 like 1		

For research use only.

GENE DESCRIPTION	BCL2-like 1
GENE ID	598
ORIGIN SPECIES	Human
ACCESSIONS	BC019307
VECTOR	pDONR223
NCBI SYMBOL	BCL2L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctcagagcaacccggagctggtgtgactttctctctacaagcttccagaagagatacagctggagtcagtttagtgatgtggaagagaacaggactgagcccccagaaggagactgaatcggagatggagaccccccagtgccatcaatggcaacccatcctggcacctggcagacagccccgcggfgaatggagccactggccacagcagcagtttgatgccgggaggtgatcccatggcagcagtaaaagcaagcgtgaggaggcagggcagcagtttgaactcgggtaccggcgggcaatcagtgacctgacatcccagctccacatcacccaggagacagcatatcagagctttgaacaggtagtgatgaactcttcgggatggggtaaactgggtcgcatgtggcctttttctccttcggcggggcactgtgcgtggaagcgtagacaaggagatgcaggtattggtgagtcgacgttgatggccacttacctgaatgaccacctagagccttgatccaggagaacggcgctgggatactttgtggaactctatgggaacaatgcagcagccgagagccgaaaggccaggaaacgttcaaccgtggttcctgacggcagcagtggtggccggtggttctgctgggctcactcttcagtcggaaatga