



Reference: PTXBC019330

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PODLX2; podocalyxin-like protein 2; podocalyxin like 2		
For research use only.			

GENE DESCRIPTION	podocalyxin-like 2
GENE ID	50512
ORIGIN SPECIES	Human
ACCESSIONS	BC019330
VECTOR	pDONR223
NCBI SYMBOL	PODXL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]