



Reference: PTXBC019341

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein transport protein Sec24A; SEC24 family, member A; SEC24 related gene family, member A; SEC24-related protein A; SEC24 homolog A, COPII coat complex component		

For research use only.

GENE DESCRIPTION	SEC24 family, member A (<i>S. cerevisiae</i>)
GENE ID	10802
ORIGIN SPECIES	Human
ACCESSIONS	BC019341
VECTOR	pDONR223
NCBI SYMBOL	SEC24A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]