

TTYH1-tweety homolog 1 (Drosophila) Gene

Reference: PTXBC019358

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein tweety homolog 1; hTTY1; tweety homolog 1; tweety family member 1		
For research use only.			

GENE DESCRIPTION	tweety homolog 1 (Drosophila)
GENE ID	57348
ORIGIN SPECIES	Human
ACCESSIONS	BC019358
VECTOR	pDONR223
NCBI SYMBOL	TTYH1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggctggcctacgtcctcctgctcctggagctgctggtctgctccttcacccctcctggcctggcgaagcagagcaagtggctggtgatcgtgatgacagtcagtgagctcctggtctcctgctcctgagctgggctccatgggctggagcagccagcgccgtgggcctcagtgac
cttctgctccaatccagaccttatgttctgaacctgacctgagggagagacagggtcagctcagacatcctgagctattatcctcctgcaaccgggcccgtctccaaccccttccaacagaggctgactctgtccagcgagctctggccaatccactccagctgctggcctggagc
gagaagctgtgctcagttcccttcagcgcagaagcctctgctgcttggaggagactctgaatgtgacagaagaaattccaccagttggtggcactgctacactgccgagcctgcacaaggactatggtgagccctgcggggcctgtgcgaagacccctggaaggcctgctct
cctgctgctcttccctgctgctcagggagcgtggccactgcccctctgagcctgcccagcctgggccccttcccacccagtgacgactacgatgacacagacgatgacgacctttcaacccctcagcaggaaatccaagcgtttgtgagtgagcagtgctctatctga