

ZMAT4-zinc finger, matrin type 4 Gene

Reference: PTXBC019598

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger matrin-type protein 4; zinc finger matrin-type 4		
For research use only.			

GENE DESCRIPTION	zinc finger, matrin type 4
GENE ID	79698
ORIGIN SPECIES	Human
ACCESSIONS	BC019598
VECTOR	pDONR223
NCBI SYMBOL	ZMAT4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagtcctccgatattgatcaggatttattcacagacagttactgcaagggtgtgcagtcgacagctgatccgaatcgagcgtgtggccactacgagagtcgaaaacatgcaagcaaatccgactgtattacatgcttcacccagggatggagggtgtcctgccaagaggctccg
gtcagaaaatggaagtgtatccgacatgggtgataagaacaagtgtgcacactctgcaacatgtcattcacttcagcgggtgtggccgattccattatcaaggcaaatccacgcaagaggttaaaactcttctaggagagaagacccattaaagaccacaggtctgaggcgcaat
tacagatgtgccatctgcagtgctccctaaactcaatagaacagtatcatgccatctgaaaggatctaaacaccagaccaacctgaagaataagtag