



Reference: PTXBC019612

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf75; NLRR-4; NLRR4; dJ1056H1.1; leucine-rich repeat neuronal protein 4; neuronal leucine-rich repeat protein 4; leucine rich repeat neuronal 4		

For research use only.

GENE DESCRIPTION	leucine rich repeat neuronal 4
GENE ID	164312
ORIGIN SPECIES	Human
ACCESSIONS	BC019612
VECTOR	pDONR223
NCBI SYMBOL	LRRN4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]