

COL20A1-collagen, type XX, alpha 1 Gene

Reference: PTXBC019637

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	bA261N11.4; collagen alpha-1(XX) chain; collagen, type XX, alpha 1; collagen-like protein; collagen type XX alpha 1 chain		

For research use only.

GENE DESCRIPTION	collagen, type XX, alpha 1
GENE ID	57642
ORIGIN SPECIES	Human
ACCESSIONS	BC019637
VECTOR	pDONR223
NCBI SYMBOL	COL20A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagaggcctggagggaactgctggcctgcctggaccctggccccaggggttcagggtcagggccaggggcactagtgagagcaggacctcaggggaccgtggggccacaggactgccaggggccaaaggggaacgaggagagaagggcgagccgca
gtccctggccacctctaccagctgtgagccaggtgtgagctgccaatcagacacacgtgtcaaaagttcgaactctccacgagaacaccagggccccatgccatcttgagcagaagctggagccgggcactgagccccctgggtccctggcaccgcagcaaggccctgg
ttcctggagaatggggcgctgtgtggccacctgagggcagaggggagcctggagctgtgtgtcagatgggcagccctgggcagcaggggctagcaccaggccctctgggagtgga