

LRRC28-leucine rich repeat containing 28 Gene

Reference: PTXBC019704

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	leucine-rich repeat-containing protein 28; leucine rich repeat containing 28		

For research use only.

GENE DESCRIPTION	leucine rich repeat containing 28
GENE ID	123355
ORIGIN SPECIES	Human
ACCESSIONS	BC019704
VECTOR	pDONR223
NCBI SYMBOL	LRRC28
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgctcgaaactttgaagacgatctgtggcaaggctagaaaagcacaagaattgttcttaaatataggaatctgcacattttccattggagtactgaaagatgagggaactgcagtacttgagagactctatatgaaaaggaaactccctgacatccttgccagaaaacctgtcag
aagcttccaaaccttgggaactgaagacaactgttagtttggcttactgtgtccaccagtgctgcagactttgacctgtgagttga