

C15orf40-chromosome 15 open reading frame 40 Gene

Reference: PTXBC019820

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0235 protein C15orf40; chromosome 15 open reading frame 40		
For research use only.			

GENE DESCRIPTION	chromosome 15 open reading frame 40
GENE ID	123207
ORIGIN SPECIES	Human
ACCESSIONS	BC019820
VECTOR	pDONR223
NCBI SYMBOL	C15orf40
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctaagaagcgtggtgcgacgaagggtaaaagccagagcaagggaaccagagagaccattctccttaggtcctgtggcagttgacctaaggatgcgtcaccatagccatccatgcaaaacctggctcacaacaaatgctgtaacagattgacagcagaggctgtaaat
gtagctattgcagcacctccatcagaggagagagctaatgctgagctctgtcggtatcttccaaggctcctagaactcaggaagagtgatgtggtttggataagggtggtaaatctcgtgaaaagggtggaagctttggcttctacaactccagaagagatctggagaaataaaaagga
agccaaaaaacataa