

ZNF511-zinc finger protein 511 Gene

Reference: PTXBC019897

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein 511		
For research use only.			

GENE DESCRIPTION	zinc finger protein 511
GENE ID	118472
ORIGIN SPECIES	Human
ACCESSIONS	BC019897
VECTOR	pDONR223
NCBI SYMBOL	ZNF511
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagttgccccgcgctgtgcgcccgcctcgtgcggggcccgccggcgagccactgcctgtagagcgggatcccgccggtggggccgccccttcgttcgttcgcgcgccccgtgcgctcccgccggagcaccagttcttcgaggatgggacgtgcagcgcacact
ctacctccaggacgtgatcatgcagggtggccgacgtgctgagaagccaggggtgcccgcgtttgctcctccaggtggccgctgctgccaggtgttcgatccctggacgactacgagcaccactaccacacgtgcacggaaatgtttgctccttttgaagcgggcttccctccgg
acacctgctggaccccacatcctgagtggtgacgattcgtcttccagatcctgtctgagggcaggacatgtatcaggtcttgtagaaggctgcacagagaagttcaagaccagcagagaccggaaggatcacatggtgaggtgacacgtgtacccccggaactccgtttgataa
gccaaagaaaagcagaagcccgccctcagcagaagcccggggacagtgagagcgggtcagaaggggagggccatggaaatctgctctgagcctgtggcagcctccctgcaccggcaggtgagagggcggtatcacagacatagaataccctctaccatctgctttgtcaggggtg
ccgctcgaggatftaaaagcaacaagaagaaaccaacaacaatgctga