



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

LPHN1-latrophilin 1 Gene

Reference: PTXBC019928

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LPHN1; C1RL1; CL1; LEC2; adhesion G protein-coupled receptor L1; C1RL-1; calcium-independent alpha-latrotoxin receptor 1; latrophilin-1; lectomedin-2		

For research use only.

GENE DESCRIPTION	latrophilin 1
GENE ID	22859
ORIGIN SPECIES	Human
ACCESSIONS	BC019928
VECTOR	pDONR223
NCBI SYMBOL	LPHN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctgctgtttccaccagtagagggtgatggcgcaagaagaaattgggaagcacatgggggttfgtgagggcacgttatggctgaggagggtgctgggaacgcgaagcgcgtctggggatggggaggggcaaaggfagagataccttcctatcctccactcccgcagggaagacgaagccagggtgccagggtctggggctggggacaccatgttctggggcggctgacctccaacaatccagaccgccgaagctcagggtcgggatgtgtctgagacactgtccggccccctctctccaggctgtcagggcacgtccccagcccccttccctaactctaggggccatagtagccccgcacccttagccccagccccctgcatfctggctttccccagagacggccccctccacctctgtctgaccttccccatctggggccccctcttgtctcacaagtgcattcgaattcacaacactctgggtcagactgtctgtcgtccatgaactctgttgggttaaggggagggggcctctgtaggaaggaactctggggcgaagggcgacaccttttgcctgctgacgt