

C19orf40-chromosome 19 open reading frame 40 Gene

Reference: PTXBC020247

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C19orf40; Fanconi anemia core complex-associated protein 24; Fanconi anemia-associated protein of 24 kDa; Fanconi anemia core complex associated protein 24		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 40
GENE ID	91442
ORIGIN SPECIES	Human
ACCESSIONS	BC020247
VECTOR	pDONR223
NCBI SYMBOL	C19orf40
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaaagaacccctgatgacgggccccgtgcacgtgcctttggggcatattgtggccaatgagaaatggcgcgggtcacagctggcgcaggagatgcaagggaataagctcatttcgaggatggctgacaccagactttatctgtcgaacagatgctgcattctttatgc
accgaagctgatttggggcaggaatggctacagaaaaggctgttgggttagaattccaataatcttaaggaaattgtagctgttgaataaacccggatgagtgaaacaatactccagccctacagaagttaactgtgctggacctggaaatgggtgctgtccagtgccagccaga
tggaagcatcctgctcgtcatccagttggttcaagagcaaaccaagagccagtaagaacctcttctcggaagaacggccctgctgtctgagccttcgctccttgaaccgtgcagcagatccaggagttgaaaagttaagctccctctcctcagaagtccaagca
tccagcaactgagtaatgcttcattggggaactggagcaggtgtcggacaagcagtgccacagcagatccatgccttcttcacgcagccaggtga