

BLOC1S2-biogenesis of lysosomal organelles complex-1, subunit 2 Gene

Reference: PTXBC020494

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BLOS2; BORCS2; CEAP; CEAP11; biogenesis of lysosome-related organelles complex 1 subunit 2; 11 kDa centrosome associated protein; BLOC-1 subunit 2; centrosomal 10 kDa protein; centrosome protein oncogene; biogenesis of lysosomal organelles complex 1 subunit		

For research use only.

GENE DESCRIPTION	biogenesis of lysosomal organelles complex-1, subunit 2
GENE ID	282991
ORIGIN SPECIES	Human
ACCESSIONS	BC020494
VECTOR	pDONR223
NCBI SYMBOL	BLOC1S2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttccaaaatggccactacctgactgggaactgacggcaccagtgaagactataagctcctggaaaatgaataaactcaccagcttgaagtatcttgaatgaaagatattgctataaacattagtaggaactaaaggacttaaccagaaatgctggactgcagccttatctg
gatcagatcaatgtcattgaagcaggtagcagctcttgagcaggcagcttacaagttggatgcataatcaaaaaaactggaagccaagtacaagaagctggagaagcgatga