

## KPNA6-karyopherin alpha 6 (importin alpha 7) Gene

Reference: PTXBC020520

|                |                                                                                                                                                          |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                      | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                     | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | IPOA7; KPNA7; importin subunit alpha-7; importin alpha 7 subunit; importin-alpha-S2; karyopherin alpha 6 (importin alpha 7); karyopherin subunit alpha 6 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | karyopherin alpha 6 (importin alpha 7)                                                           |
| <b>GENE ID</b>          | 23633                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC020520                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | KPNA6                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggagaccatggcgagccagggaagacaattatcgaatgaagagctataagaacaattcttaaccctgaagaaatgagacgaagaagagaggagggcattcagctccggaagcagaagcgagagcaacaactttttaacggagaatgtggagctgattaatgaagaa  
gctgccatgttcgatagctcttctcatggactcttatgtgagctctaccactggggagagtgatcacaaagagagatgggtggagatgctcttttctgatgattctgacctgcagtttagcaaccacacagaaatccggaaactgctctccaaagagccttagtctccaatagatgaagttatcaaca  
ctccaaagatgggtgagcgggttcgtggagttctgaagaggaatgagaattgtacattacagtttgaagctgctgggctctaaacgaatattgctctggaacctctcagcagaccaaattgtcattgaagcaggggctgtcccatttttatagagctgcttaattcagactttgaggatgttcag  
gaacagggcagctctgggcaactgggaacatagctggagatagctctgttttccgagattacgtcttgaactgttccatcttaactcttggtaacactcttaccagtgccacacgactgacgatgacacggaatgcagtgctgggcccgtgtcaaatctctgccagggaacaccccccag  
agtttgcagggtctctctgttttgcctgtactgtctcgcctactcttcagcagcagctcggactgtgctggcagatgctgtggcccctttctatctgtctgtatggcccaatgagaagatccaggcagtcatagactccggagctgctggagattggtagagctgctgatgcacaatgatta  
caaatgtgcttctcctgagagccgtgggttaacatgctcactgggatgacatccagaccaggtcattcttaactgttcagccctacctgcttctccactgtttgagcagtgcccaaggagtcattccggaaggaaagcttgcctggactattcaaatattactgtgccaacagggtc  
aaatacaggctgttatagatgcaaatatcttccctgtgtgtgatcgaatcttcagaaagcagagtttctgtaacaggaaagaggcagcctgggccaatcaccaatgccacatcaggaggaacccctgagcagatcaggttactgtgctcactgggctgcatcaaacccctatgtgactgtctga  
ctgtaatgattcgaagattgtcgaagtggccctcaatggactggagaacatctgcggcttgagagcaagagggaagcagctggtcaggggtcaatcttattgtggcctcatagagggaagcctatgcttggataaaatgagtttctccagagccacgagaaccaggagatcta  
ccagaaggccttcgactcattgagcactatttgggttagaagacgatgatagcagctggctcccaagtcgatgaacgcaacagcagtttcatcttcagcagcctgaggccccatggaggcttcagctataa