



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C12orf23-chromosome 12 open reading frame 23 Gene

Reference: PTXBC020522

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

UPF0444 transmembrane protein C12orf23; C12orf23; transmembrane protein 263

For research use only.

GENE DESCRIPTION	chromosome 12 open reading frame 23
GENE ID	90488
ORIGIN SPECIES	Human
ACCESSIONS	BC020522
VECTOR	pDONR223
NCBI SYMBOL	C12orf23
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]