

NIP30-NEFA-interacting nuclear protein NIP30 Gene

Reference: PTXBC020536

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NEFA-interacting nuclear protein NIP30; NIP30; C16orf94; CDA018; CDA10; protein FAM192A; family with sequence similarity 192 member A		

For research use only.

GENE DESCRIPTION	NEFA-interacting nuclear protein NIP30
GENE ID	80011
ORIGIN SPECIES	Human
ACCESSIONS	BC020536
VECTOR	pDONR223
NCBI SYMBOL	NIP30
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatggaggggatgatgtaacctattatcaaaaagaggtttgtgtctgaggcagaactagatgaacgcgcgcaaaaggaggcaagaagaatgggagaaagttcgaaacctgaagatccagaagaatgtccagaggaggtttatgacctcgatctctatatgaaaggctacaggaa
cagaaggacaggaagcagcaggagtagcaggaacagttcaaatcaaaaacatgtaagaggttagatgaagatgagaccaactctctgaggtttctcgacagcagggaactaatagaaagcaacgaagagaagaagaaactgaagaactgaaggatacagaataacctca
agaaggttggaatttctcaagagaacaagaaggaaagtggaagaaactgactgtgagcctatagaacaaagaacaagttctccaggcggaagctgttggcaggagctgtgaagcataagagctcagagagtggaacaggtgtgaaagactgaaccggaccctgagccagatg
acaagaatcaagagccctcatctgcaagtcctcggaaacacctccctgagtgcccccctccatccactgcccctctgctgcagtatgtatcgcatctccaggccctgggtgctactctgggagcagcgaactccaggtccagctcagacagcgaaggcaccatcaatgccaccgaa
agattgtctctccatctccgaaccaacacctctctgagggcccctag