

SNX12-sorting nexin 12 Gene

Reference: PTXBC020559

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	sorting nexin-12; sorting nexin 12		

For research use only.

GENE DESCRIPTION	sorting nexin 12
GENE ID	29934
ORIGIN SPECIES	Human
ACCESSIONS	BC020559
VECTOR	pDONR223
NCBI SYMBOL	SNX12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggacacgagtagctgataccggcgcttaactcgaagccgcaggacctgaccgacgttacggcgcccaagtaacttctggagatcgacatcttaatcctcagacgggtggcgctgggacgcgcgcttcaccacctatgaggttcgcatgcggacaaacctacatc
ttcaagctaaaggagtcctgctacggcgctacagtgactttgagtgctgaaatgagctggagagatagcaagattgtatgaccacctgcctgggaagccctgaagcggcagctcccttccgaggagatgaaggatctttgaggagctttcatcgaagaaggagg
caggcctcgagcagtttatacaaaattgctggcaccactggtcagaatgaacgctgcctacacatgttctgcaagaggaggcaattgacaggaaactacgtccggggaagagcttctgtgtcttcccaggctggagtgcatggcatga