



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

HAMP-hepcidin antimicrobial peptide Gene

Reference: PTXBC020612

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

HEPC; HFE2B; LEAP1; PLTR; hepcidin; liver-expressed antimicrobial peptide 1; hepcidin antimicrobial peptide

For research use only.

GENE DESCRIPTION	hepcidin antimicrobial peptide
GENE ID	57817
ORIGIN SPECIES	Human
ACCESSIONS	BC020612
VECTOR	pDONR223
NCBI SYMBOL	HAMP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcactgagctcccgacatctgggcgcttgcctctgtcctctctctctcgcgcagcctgaccagtggctctgtttcccacaacagacgggacaaccttgacagctgcaacccccagcacagagctggaggccagggccagctggatgcccatgtlccagaggcgaaggaggcgag
aacaccacttccccctgcatttctctgcgcgctgctgtcatgatcaaagtgtggatgtctgccaagactgag