

C8orf4-chromosome 8 open reading frame 4 Gene

Reference: PTXBC020623

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein C8orf4; TC-1; TC1; human thyroid cancer 1; thyroid cancer protein 1; chromosome 8 open reading frame 4		

For research use only.

GENE DESCRIPTION	chromosome 8 open reading frame 4
GENE ID	56892
ORIGIN SPECIES	Human
ACCESSIONS	BC020623
VECTOR	pDONR223
NCBI SYMBOL	C8orf4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagcaagcggaagccaccaagccatcatgtccacgtcgtacgagtcagcccatcatcctggtaccattcgacacagcctctgtaagaaagccgtgggcaacatcttgaacacagacagcaagaatcactagaaggtctctcagaactctggagacaagaagcagagagagagccaagatctttgcatagatcaagatgtggaggagaaaacgctgcccctgatggccttgaagaagaggacaaaagacaagctttccagttctgaaactgcggaaatattccatcaaaagttcactga