

CCDC59-coiled-coil domain containing 59 Gene

Reference: PTXBC020647

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BR22; HSPC128; TAP26; thyroid transcription factor 1-associated protein 26; TTF-1-associated protein 26; TTF-1-associated protein BR2; coiled-coil domain-containing protein 59; coiled-coil domain containing 59		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 59
GENE ID	29080
ORIGIN SPECIES	Human
ACCESSIONS	BC020647
VECTOR	pDONR223
NCBI SYMBOL	CCDC59
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccggtgagcggtccggaagtggcgccctggtgtattgagcgcggtggaaggggtttccactgtcgggtacaggaataagaatgtgagacagaagacatggcgccctaaccacccgcaagcctctgtggggagcgttcgcgagggacaaggcttgcctttcgaaga
aaactgaaaatacagcaaaagttaacaagaattgctacggaaggaaaagaagctcaaacgtcactggaatctcaatfcacagatcgataccagataatctgaaacatctctatttagctgaagaggaaagacataggaagcaagcaagaaaagtcgaccatcttctcagaacaagtcca
ccagccgttgcctgaagaacagtgtagcattgacgagccttatttgaagatcagtgtagctttgaccagcctcagccagaagaacaatgtattaaaacagtaaacctcttacaattccaagaaaaataaaaaaagaaacatcaaatcaaaagcacaagaagaatgaacagatacagcc
aaacgtgctgctaagaacaagaattcgagaggagaaaacaggagagagaagaagcccaaggcagtaaaaaagaagaatggaagtgtttaaatactgaacaaaagactaaaaagggccaaccaactgaatgtacaatggagtccttctcaaaaaatacaagaaaat
gttaa