

HIGD1B-HIG1 domain family, member 1B Gene

Reference: PTXBC020667

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CLST11240; CLST11240-15; HIG1 domain family member 1B; HIG1 hypoxia inducible domain family member 1B		

For research use only.

GENE DESCRIPTION	HIG1 domain family, member 1B
GENE ID	51751
ORIGIN SPECIES	Human
ACCESSIONS	BC020667
VECTOR	pDONR223
NCBI SYMBOL	HIGD1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgctaacagacgctgggtaccacctgacgacgaagactgtgtgtctgagaagctcctgaggaagactcgggaatccactgtgtcctataggcttaggaggtgctgtgtgtgtagcagcatacaggattaccggctgaggtctcgtggtccaccaagatgtccatacact
gattcacaccgagtgagcagcgaggcctgtgcagtggtgcaatcatgctaggtgctgtgtacacaatgtacagcgattacgtcaagaggatggcacaggatgctggagagaagtag