

SLPI-secretory leukocyte peptidase inhibitor Gene

Reference: PTXBC020708

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ALK1; ALP; BLPI; HUSI; HUSI-I; MPI; WAP4; WFDC4; antileukoproteinase; HUSI-1; WAP four-disulfide core domain protein 4; mucus proteinase inhibitor; protease inhibitor WAP4; secretory leukocyte protease inhibitor (antileukoproteinase); seminal proteinase i		

For research use only.

GENE DESCRIPTION	secretory leukocyte peptidase inhibitor
GENE ID	6590
ORIGIN SPECIES	Human
ACCESSIONS	BC020708
VECTOR	pDONR223
NCBI SYMBOL	SLPI
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagtcacagcgccctctccctctctggtgctgctgcccctgggaactctgacaccttgggtgtggaaggctctggaagtcctcaaaagctggagctgtctctcctaaagaaatctgccagtgcccttagatacaagaacctgagtgccagagtgactggcagtgccaggggaagaa
gagatgtgtcctgacactgtggcatcaaatgctggatcctgttgacaccccaaaccaagaaggaggaagcctgggaagtgccagtgacttatggccaatgtttgafgcttaaccccccaatttctgtgagatggatggccagtgcaagcgtgacttgaagtgttgcatggcatgtgt
gggaatcctgcgttccctgtgaaagcttga