

CRP-C-reactive protein, pentraxin-related Gene

Reference: PTXBC020766

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PTX1; C-reactive protein; C-reactive protein, pentraxin-related; pentraxin 1		

For research use only.

GENE DESCRIPTION	C-reactive protein, pentraxin-related
GENE ID	1401
ORIGIN SPECIES	Human
ACCESSIONS	BC020766
VECTOR	pDONR223
NCBI SYMBOL	CRP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaagctgtgtgttcttggtcttgaccagcctctctcatgcttttggccagacagacatgctgagggaaggcttttgttttccaaagagtcggatacttctatgtatccctcaagcaccgttaacgaagcctctcaagccttcaactgtgtgcctccacttctacacggaactgtcctcg
accctgtgtcctaatgtcctgaactggcgggcactgaagtatgaagtgcaggcgaagtgtcacaaacccagctgtgcccctga