

SAPS3-SAPS domain family, member 3 Gene

Reference: PTXBC020953

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SAPS3; C11orf23; PP6R3; SAP190; SAPL; SAPLa; serine/threonine-protein phosphatase 6 regulatory subunit 3; SAPS domain family, member 3; sporulation-induced transcript 4-associated protein SAPL; protein phosphatase 6 regulatory subunit 3		

For research use only.

GENE DESCRIPTION	SAPS domain family, member 3
GENE ID	55291
ORIGIN SPECIES	Human
ACCESSIONS	BC020953
VECTOR	pDONR223
NCBI SYMBOL	SAPS3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggtggacctgagtgaaacccaactggtcagtaactttgatgtcccaatggaacaacccacgggtgctcattggattctgtgggatctgatgtctggagcacagaggagccgatgccaaactaaagacgggctgggctcttttcagagttcacgtcttcctgagcaca
gattctttaaggagtaattctccagtggaatggaaccagcactgaacccatggaccctctgactccagtgccgtgcccctggcagtgccagcagaagcggcagtgctggtccatggaagccagctctgacggagaggagatgcagaaagtacagacaaggttaactgagaca
gtgatgaatggcgcatgaaggaaacgctcagcctcactgtagatgccaagacagagactgctgtctcaaaagtggaggaaggaaactgtctacctctcaagatgctgtgttaagacgcagaggagtgctcccgagactgcagaggcgaagtgcggcgccagcctccagc
agcagctccgagcagaggactggccaaccaagcgcaccaggtgacacttcagtgaatggcctgtatga