

SEMA4G-sema domain, immunoglobulin domain (Ig), transmembrane domain (TM) and short cytoplasmic domain, (semaphorin) 4G Gene

Reference: PTXBC020960

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	semaphorin-4G; sema domain, immunoglobulin domain (Ig), transmembrane domain (TM) and short cytoplasmic domain, (semaphorin) 4G; semaphorin 4G		

For research use only.

GENE DESCRIPTION	sema domain, immunoglobulin domain (Ig), transmembrane domain (TM) and short cytoplasmic domain, (semaphorin) 4G
GENE ID	57715
ORIGIN SPECIES	Human
ACCESSIONS	BC020960
VECTOR	pDONR223
NCBI SYMBOL	SEMA4G
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggctcaatgtcaccacctctgcatggccctgtgtgctggtgctgtgaaaccagacaagacctctgccagccacctaagccctgcgtacattcacatgcacacatggaagaatgtttatggctgggctgcagtccccaccctcacctctctctggtgcattctgtttcatccctgctctggacttgggtaccctcccaattgccacatcctatctgtctctctccagccccatgtgtgacctctttgtcaagagcttgggaacggccagcctggggaggtaagactgcatcactcccctctctccctctctgtgtggcccttgaatcagcctccccactctctgtgctattctcaagagtatga